



## HIGH PRECISION SYRINGE PUMPS

10905 Cash Road, Stafford, TX 77477  
Phone: 281.277.5499 / Fax: 281.277.0045  
sales@chemyx.com / chemyx.com

# Fusion 4000 Updated Serial Commands for Full Computer Control

help

-----  
help            -Show Help Information  
-----

## Movement Commands

-----  
[x] start [y]            - Pump [x] [y](0 Basic/1 Program) Run  
[x] stop                - Pump [x] Stop  
[x] pause               - Pump [x] Pause  
start [y]               - All Pump [y](0 Basic/1 Program) Run  
stop                    - All Pump Stop  
pause                   - All Pump Pause  
restart                  - Pump Restart  
-----

## Setting Commands

-----  
[x] set units [xxx]      -Pump [x] Set Units(0,1,2,3 or ml/min,ul/hr)  
[x] set diameter [x.x]   -Pump [x] Set Diameter  
[x] set rate [x.x,x.x,...] -Pump [x] Set Rate  
[x] set volume [x.x,x.x,...] -Pump [x] Set Volume  
[x] set delay [xxx]      -Pump [x] Set Delay  
[x] set time [x.x,x.x,...] -Pump [x] Set Time(will calculate rate)  
-----

## Status Commands

-----  
[x] read limit parameter [y] -Pump [x] Return [y](0 Basic/1 program) Max/Min Volume & Rate  
[x] dispensed volume      -Pump [x] Return Volume Dispensed  
[x] elapsed time           -Pump [x] Return Time Elapsed  
[x] pump status            -Pump [x] Return System Status  
-----

## Additional Optional Commands

-----  
[x] hexw2 [units] [mode] [diameter] [volume] [rate] [delay]      -Sends All Settings  
[x] hexw2 [units] [mode] [diameter] [volume] [rate] [delay] start   -Sends All Settings and Starts  
-----

The red part will be the updated commands.



## HIGH PRECISION SYRINGE PUMPS

10905 Cash Road, Stafford, TX 77477  
Phone: 281.277.5499 / Fax: 281.277.0045  
sales@chemyx.com / chemyx.com

- [x] can be 1, 2, or 3
  - When [x] = 1, that controls Pump Channel 1
  - When [x] = 2, that controls Pump Channel 2
  - When [x] = 3, that controls Pump Channel 1,2 Cycle Mode
- [y] can be 0 or 1
  - When [y] = 0, that sets the pump in Basic Mode
  - When [y] = 1, that sets the pump in Programmable Mode

For example,

- 1 start 1\r = start Pump 1 in Programmable Mode
- 2 set rate 5 = sets Pump 2 to infuse at a rate of 5 unit/time
- 1 set volume 1.0,-2.0,3.0\r = sets Pump 1 to Programmable Mode to:
  - First step volume: infuse 1.0
  - Second step volume: withdraw 2.0
  - Third step volume: infuse 3.0
- 1 set volume 2.0 loop 1, -1.0 loop 2\r = sets Pump 1 to loop in Programmable Mode
  - First step volume: infuse 2.0; loop once
  - Second step volume: withdraw 1.0; loop twice

To set a step to withdrawal mode, add a negative value before the volume!