

CHEMYX

STAINLESS STEEL SYRINGES SPECIFICATIONS

SCALE VOLUME (ml)	INNER ID (mm)	MAX. PRESSURE (bar)	WEIGHT (lb)	Φ x LENGTH (mm)
6	9.53	150	0.36	13.5 x 167
10	12.60	110	0.5	17 x 167
20	19.13	50	1.0	23 x 167
50	28.60	20	1.34	33 x 167
100	34.90	10	2.0	39 x 189
200	43.96	5	3.30	49 x 206
300	54.00	5	4.19	59 x 206

