


OEM Module - Data Sheet



The OEM syringe pump is available in both Full Case Model 30060 and Closed Case Model 30050

TECHNICAL SPECIFICATIONS

MODE	Infuse Only
MINIMUM FLOW RATE	0.0001ul/min (with 0.5ul syringe)
MAXIMUM FLOW RATE	140mL/min (with 50mL syringe)
LINEAR FORCE	35lb (15.9kg)
SYRINGE SIZE (MIN/MAX)	0.5mL to 60mL
STEP RESOLUTION	0.098 μ m/step
OPERATIONAL TEMPERATURE	4°C to 40°C (40°F to 104°F)
STORAGE TEMPERATURE	-10°C to 70°C (14°F to 158°F)
ACCURACY	\pm <0.35%
REPRODUCIBILITY	\pm <0.05%
CONNECTIVITY	USB-B, RS232 and TTL Ports
POWER	110/220 VAC

SELECT APPLICATIONS

- Microfluidics
- Animal Infusion
- Neuroscience
- Mass Spectrometry Infusion
- Electrospinning
- Reagent Dispensing
- Microdialysis
- Stem Cell Injections
- Reactor Feeding

GET MORE - For more tips, a catalog of support articles, and step-by-step tutorial videos, visit www.chemyx.com/support.